


Cortex Gas Analyser User Guide




Initial Setup

- Power on Cortex 15 minutes prior to calibrating to warm up sufficiently
- Equipment required:
 - Face mask
 - Head net
 - Heart rate monitor
 - Calibration gas canister
 - Mouthpiece turbine


Metasoft Setup

- Load Metasoft 3 program, icon on desktop 
- Select calibration tab



Pressure Setup

- Select Pressure Tab
- Enter Room Pressure from barometer
- Select Start
- Select Transfer
- Close once complete



Gas Calibration

- Select Gas tab

MetaSoft 3.9

Medical Diagnostics

Database

Calibration

Pressure

Gas

Volume

Check

Edit

Print

O2 Sensor Exchange

Configuration

Local Operator: Administrator Device Configuration: cortex

Gas Sensor Calibration

Settings

Gas 1
Ambient Air (O2: 20.93 / CO2: 0.03 Vol%)

Gas 2
cal gas (O2: 17.05 / CO2: 4.98 Vol%)

Edit Span Gases...

Calibration Measurement

Time

Start Gas 1 Start Gas 2

Results

	O2	CO2
Factor	<input type="text"/>	<input type="text"/>
Offset	<input type="text"/>	<input type="text"/>

Transfer

Help Close

Gas calibration

- Place gas sensor in open air
- Select Start Gas 1, dots will appear on the two lines and wait until calibration is complete



Gas Sensor

The screenshot shows the MetaSoft 3.9 software interface. The main window has a menu bar (MetaSoft, Help) and a sidebar with various diagnostic tools. The 'Gas Sensor Calibration' dialog box is open, displaying the following information:

Settings

Gas 1: Ambient Air (O2: 20.93 / CO2: 0.03 Vol%)
Gas 2: cal gas (O2: 17.05 / CO2: 4.98 Vol%)

Calibration Measurement

Time 0:00:40

The graph shows two horizontal lines: a blue line for O2 and a pink line for CO2. Three data points are plotted on each line.

Results


	O2	CO2
Factor	<input type="text"/>	<input type="text"/>
Offset	<input type="text"/>	<input type="text"/>

Buttons: Stop, Start Gas 2, Transfer, Help, Close.

System tray: Local, Operator: Administrator, Device Configuration: cortex


Gas Calibration

- Connect gas sensor to gas bottle
- Select Start Gas 2



Gas Calibration

- Once complete select transfer and close



Turbine and Sensor Connection

- To connect both parts, twist the bottom of the turbine until the small ridge at the bottom is almost totally obscured
- Then slide in the sensor and twist the bottom of the turbine once again to lock it in place. Note the groove in the sensor needs to be facing the bottom of the turbine



Locking groove

Sensor



Twist this part

Turbine

Volume Calibration

- Connect Hans-Rudolph device to gas turbine and sensor




Volume Calibration

- Insert plug into gas sensor before calibrating volume



Volume Calibration

- Ensure that the inspiration (pull lever out) and expiration (push lever in) are kept between the green boundaries



- 5 good strokes need to be registered and need to be as close to the green line for the most accurate readings
- Once completed select Transfer and close


Connecting Mask and Turbine

- Remove the plug from the gas sensor
- Connect turbine to a face mask and head net
- Plug in gas sensor
- Ensure that the cables come out of the top of the mask to avoid blockages



Adding a Test Subject

- Select the Medical Diagnostics tab
- Select CPX Testing



Adding a Test Subject

- Select Prepare
- Select New Test Subject

The screenshot shows the MetaSoft 3.9 software interface. The 'Prepare' menu is open, and 'New Test Subject...' is selected. The interface includes a menu bar, a toolbar, a 'Test Subject:' field, a 'PULMONARY MONITOR' graph, a 'WORKLOAD DEVICE CONTROL' panel, and a data display table.

PULMONARY MONITOR

WORKLOAD DEVICE CONTROL

No Parameter	V'O2/kg [ml/min/kg]	V'E [l/min]	BF [1/min]
	0	0.0	0
V'E/V'CO2 []	V'CO2 [l/min]	V'O2 [l/min]	RER []
0.0	0.0	0.0	0.0

Local | Operator: Administrator | Device Configuration: cortex

Adding a Test Subject

- Enter Test Subjects details
 - Name
 - Sex
 - Date of birth
 - Height (cm)
 - Weight (kg)
 - Physical Activity
 - Cortex mask size
- Select save and close

MetaSoft 3.9

Medical Diagnostics Measurement: CPX Test

CARDIOPULMONARY MONITOR WORKLOAD DEVICE CONTROL

New Test Subject

Standard | Medical Baseline Information | Details/Notes | Address

Identification Number: [] Initials: [] Height [cm]*: 184

Title: [] Sex: male female Weight [kg]*: 70

Name*: Armstrong Date of Birth [dd/MM/yyyy]*: 01/01/1970 Physical Activity*: >6 hour/week

First Name*: Lance Spirometric Standard: None Mask*: CORTEX Adult 5 (medium)

Name Addition: [] Hematocrit (Hct): 0 % (Volume [ml]: 47)

Target Value Parameter: [] Teams*: Team

Target Value: 0.00 %

* Mandatory Field


Print New Save Help Close

0.0 0.0 0.0 0.0

Local Operator: Administrator Device Configuration: cortex

Adding a Test Subject

- When the test is ready to begin, select the green play button to start measurement
- An ambient air measurement will take place. This can be skipped by clicking close and then ignore in the following pop up box that appears




Extras

- On the navigation pane at the top, you can add markers which will highlight a certain point in time that may be of interest (e.g. the beginning of a test)
- The screens tab at the top can change the appearance of the data collection screen (field being the most common view – as seen in lab classes)


Saving Data to Excel

- Once a test is completed a series of boxes will appear. Close them all.
- Then select Test
- Select Export
- Select Spreadsheet interface



Saving Data to Excel

- Select the parameters you require
- Select the > to move the parameters selected in the left box across or >> to transfer all the parameters
- Excel will then open with all the data from that test
- Ensure that you click save as and change the save type to excel 97-2003 so it is compatible with your computer
- Email to yourself as USB sticks can carry viruses which can be transferred onto University PC's and cause them to run slowly



Cleaning up afterwards!

- Once you have finished testing fill a bowl with warm water and add a small amount of Neutracon detergent (10ml is enough so no need to pour a large amount in)
- Place everything in this bowl – Heart rate monitors, hair nets, face masks, turbines
- NOTE: Do NOT rinse turbines under a tap as this may break them, but flushing them in the bowl is fine
- Once washed, empty the bowl or use another bowl to disinfect the turbines and face masks ONLY, head nets and heart rate straps will lose their elasticity if left too long. Fill the bowl with cold water and add a small amount of Milton disinfectant. Then leave for 15 minutes as this is the time it takes to fully disinfect the items.
- Shake off as much water as possible as this will reduce the time it takes to dry off in the drying cabinet. NOTE: do not put anything on the very bottom and do not change the temperature as this may cause things to melt.
- Spray the handlebars/treadmill with Virkon solution (pink spray bottles) and wipe with tissue to ensure that any blood/sweat is cleaned off for the next person to use it.